EmbryoFreeze™ EmbryoThaw™

EN

Medium for freezing and thawing human embryos

Doc. reference: FP09 I10 R01 C.3
Update: 01.10.2012

INTENDED USE

EmbryoFreeze and EmbryoThaw (EmbryoFreeze/Thaw) are a set of ready-to-use media for freezing and thawing human embryos between 2PN and 4 cell stage.

COMPOSITION

EmbryoFreeze/Thaw are propandiol based cryopreservation media which also contain PBS, sucrose and human serum albumin (14-15g/Liter). EmbryoFreeze/Thaw do not contain antibiotics.

MATERIAL INCLUDED WITH THE KIT

EMF30_P* (30 procedures)

- » 3 vials of EmbryoFreeze medium (10ml)
- » 3 vials of EmbryoThaw, thawing 1 (10ml)
- » 3 vials of EmbryoThaw, thawing 2 (10ml)
- » 3 vials of EmbryoThaw, thawing 3 (10ml)

EMF40_T_P (40 procedures)

- » 4 vials of EmbryoThaw thawing 1 (10ml)
- » 4 vials of EmbryoThaw, thawing 2 (10ml)
- » 4 vials of EmbryoThaw, thawing 3 (10ml)

EMF120_F_P

- » 12 vials of EmbryoFreeze medium (10ml)

* The media should be used in the order displayed per product code.

MATERIAL NOT INCLUDED WITH THE KIT

- » sterile pipettes
- » well dishes
- » freezing straws
- » automatic freezer
- » freezing tank with liquid nitrogen (LN2)
- » waterbath (able to hold 37°C)
- » LAF bench (ISO Class 5)
- » microscope
- » labtimer

EMBRYOFREEZE/THAW AND EMBRYOCULTURE

EmbryoFreeze an EmbryoThaw can be used in combination with FertiCult media (FertiCult Flushing medium, FertiCult IVF medium and FertiCult G3 medium) before freezing and after thawing.

PRODUCT SPECIFICATIONS

- » Chemical composition
- » pH between 7,20-7,40
- » Osmolality EmbryoThaw 3: 530-560 mOsm/kg
- » Sterility: SAL 10⁻³
- » Endotoxins: < 0,25 EU/ml
- » Mouse Embryo Assay (blastocysts after 96h) ≥ 80%
- » Use of Ph Eur or USP grade raw materials (if applicable)
- » The certificate of analysis and MSDS are available upon request.

PRE-USE CHECKS

- » Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- » Do not use the product if seal of the container is opened or defect when the product is delivered.

STORAGE INSTRUCTIONS

- » Store between 2-8°C.
- » Do not freeze before use.
- » Keep away from sunlight.
- » The products can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- » Do not use after expiry date.

WARNINGS AND PRECAUTIONS

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination.

Only for the intended use.

METHOD

Ensure all media are well mixed before use.

Freeezing

1. Using a sterile pipette place 1ml of EmbryoFreeze freez-ing medium in a centre well dish (at room temperature)
2. Add the embryos (2PN stage) to the freezing medium and allow them to settle for about 30 seconds. Caution: because of density differences the embryos tend to float upwards and shrink like raisins.
3. Load the embryos in straws leaving about 1/5 air in the straw.
4. Seal the straws and label with name, date and number of embryos.
5. Start freezing program within 5 to 10 minutes. Below is an example of a freeezing protocol.

	Temperature range	Freeezing rate	Time
Step 1	Room temperature to +4°C	-10°C/min	2 min
Step 2	+4°C to -6°C	-2°C/min	5 min
Step 3	-6°C (autoseeding)	0°C/min	10 min
Step 4	-6°C to -30°C	-0.3°C/min	80 min
Step 5	-30°C to -196°C	-199°C/min	1 min

Thawing

1. Place in a 4-well culture dish 1ml of each EmbryoThaw thawing solution (1, 2 & 3). This leaves 1 well empty to retrieve the frozen/thawed embryos.
2. Prepare a water bath at 37°C to thaw the straws. Prepare a 1ml tuberculin-syringe by filling it with 0.8ml of air first followed by 0.2ml of thawing 1.
3. Remove the straws from liquid nitrogen and leave at room temperature for about 5 seconds.
4. Submerge the straw in the waterbath at 37°C for another 5 seconds (ensure no frozen part remains in the straw).
5. Empty the straw by opening both sides of the straw (above the empty well) and blowing the contents of the syringe through the straw.
6. Retrieve the embryos using a microscope and place them in thawing solution
7. Transfer the embryos to thawing solution 2 after 3 to 5 minutes.
8. After another 3 to 5 minutes the embryos are transferred to thawing solution 3. Leave for a further 3 to 5 minutes before proceeding.
9. At the final stage the embryos are transferred in IVF culture medium for washing and further culture.

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