ID: FP09 I11 SSP R01 B.2, Update: 15/11/2016
SpermFreeze SSP is sterilized by sterile filtration

GENERAL INFORMATION AND INTENDED USE

SpermFreeze SSP is a ready-to-use cryopreservation medium for freezing human sperm cells at ultra low temperatures (-196°C), in liquid nitrogen.

COMPOSITION

SpermFreeze SSP is a HEPES buffered freezing medium for use with human sperm^{1,2}. SpermFreeze SSP contains 0.4 % human serum albumin and 26.76% glycerol to protect the sperm from damage due to the freezing procedure³.

Note: The concentration of glycerol in SpermFreeze SSP is significantly higher than in SpermFreeze, it is important to use the correct dilution factors when using SpermFreeze versus SpermFreeze SSP.

MATERIAL INCLUDED WITH THE KIT

Product code	Product description
SSP001	5x 1 ml of SpermFreeze SSP (glass bottle)

MATERIAL NOT INCLUDED WITH THE KIT

- Sperm freezing straws (e.g. CBS high security sperm straws)
- Sterile pipettes
- Freezing tank with liquid nitrogen
- LAF bench (ISO Class 5)

PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20 - 7.90 (release criteria: 7.20-7.60)
- Sterility: sterile (SAL 10⁻³)
- Endotoxin: < 0.25EU/ml
- Sperm survival test: ≥ 80% survival after 4 hours exposure of untreated semen to the test medium
- Not MEA tested
- Use of Ph Eur or USP grade products if available
- A certificate of analysis is available upon request.

PRE-USE CHECKS

- Do not use the product if it becomes, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defect when the product is delivered

STORAGE INSTRUCTIONS

- Store SpermFreeze SSP between 2-8°C.
- Do not freeze before use.
- Keep away from (sun)light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and products are stored at 2-8°C
- Do not use after expiry date.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).

SPERMFREEZE SSP AND SPERM PREPARATION

SpermFreeze SSP can be used in combination with Sil-Select Plus™^{4,5}.

Before freezing

In case of very low sperm concentrations it is advisable to concentrate the sperm before freezing. This may increase the sperm quality after thawing and will reduce the number of straws to be frozen.

After thawing

If necessary, use sperm preparation techniques after thawing to eliminate dead sperm cells and debris. Dilute the concentrated sperm in washing medium or any other medium you would normally use.

METHOD

Ensure all media are well mixed before use.

Freezing

1. Allow the semen to liquefy at room temperature for 30 minutes.
2. Mix 3 parts sperm with 1 part SpermFreeze SSP. Add the medium in drops while gently swirling.
Caution: To avoid cold-shock, make sure SpermFreeze SSP is at room temperature.
3. Leave the mixture for 10 minutes at room temperature for equilibration.
4. Suck the sample/medium mixture into the freezing straws, leaving approximately 1.5cm of air at the end of the straw.
5. Seal the straws.
6. Dry off individually with a linen free cloth.
7. Shake to move the air-bubble to the centre of the straw.
8. Freeze vertically for 15 minutes, just above the level of the liquid nitrogen.
9. Store in liquid nitrogen.

Thawing

1. Remove as many straws as required from the liquid nitrogen.
2. Place the straws in tap water for 5 minutes.
3. Cut off the end of the straw, place the open end inside a container (e.g. a test tube) and tap the straw against the side of the container to allow complete evacuation of the mixture.
4. Dilute the concentrated sperm in a suitable insemination medium (at least 3 ml per 0.5 ml semen) and mix thoroughly
5. Centrifuge during 15 minutes at 300-350g
6. Resuspend pellet in a suitable insemination medium (e.g. FertiCult Flushing medium)

WARNINGS AND PRECAUTIONS

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be completely excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.

The above media do not contain antibiotics, always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid contamination.

BIBLIOGRAPHY

- 1 Mahadevan M, Trounson AD. Effect of cryoprotective media and dilution methods on the preservation of human spermatozoa. *Andrologia*, 1983; 15: 355-66.
- 2 Mahadevan M, Trounson AD, Leeton JF. Successful use of human semen cryobanking for in vitro fertilization, *Fertil Steril*, 1983; 15: 355-66.
- 3 Brotherton J. Cryopreservation of human semen. *Archives of Andrology*, 1990; 25: 181-95.
- 4 Kobayashi T, Kaneko S, Hara I, Park YJ, et al. Concentrating human sperm before cryopreservation. *Andrologia*, 1991; 23: 25-8.
- 5 Graczykowski JW, Siegel MS. Influence of sperm processing on the fertilizing capacity and recovery of motile sperm from thawed human semen. *Archives of Andrology*, 1991; 26: 155-61.

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- 6 Mahadevan M, Trounson AD. Effect of cryoprotective media and dilution methods on the preservation of human spermatozoa. *Andrologia*, 1983; 15: 355-66.
- 7 Mahadevan M, Trounson AD, Leeton JF. Successful use of human semen cryobanking for in vitro fertilization, *Fertil Steril*, 1983; 15: 355-66.
- 8 Brotherton J. Cryopreservation of human semen. *Archives of Andrology*, 1990; 25: 181-95.
- 9 Kobayashi T, Kaneko S, Hara I, Park YJ, et al. Concentrating human sperm before cryopreservation. *Andrologia*, 1991; 23: 25-8.
- 10 Graczykowski JW, Siegel MS. Influence of sperm processing on the fertilizing capacity and recovery of motile sperm from thawed human semen. *Archives of Andrology*, 1991; 26: 155-61.

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RIFERIMENTI

Mahadevan M, Trounson AD. Effect of cryoprotective media and dilution methods on the preservation of human spermatozoa. *Andrologia*, 1983; 15: 355-66.

Mahadevan M, Trounson AD, Leeton JF. Successful use of human semen cryobanking for in vitro fertilization, *Fertil Steril*, 1983; 15: 355-66.

Brotherton J. Cryopreservation of human semen. *Archives of Andrology*, 1990; 25: 181-95.

Kobayashi T, Kaneko S, Hara I, Park YJ, et al. Concentrating human sperm before cryopreservation. *Andrologia*, 1991; 23: 25-8.

Graczykowski JW, Siegel MS. Influence of sperm processing on the fertilizing capacity and recovery of motile sperm from thawed human semen. *Archives of Andrology*, 1991; 26: 155-61.

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FertiPro N.V.
Industriepark Noord 32
8730 Beernem – Belgio
Tel. +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
<http://www.fertipro.com> – E-mail: info@fertipro.com

0344