## VitriFreeze™ VitriThaw™

(EN)

Medium for vitrification and thawing of human embryos

**STERILE** **A**

Sterilised by sterile filtration.
Doc. reference: FP09 I46 R01 B.5
Update: 29.01.2018

**INTENDED USE**

VitriFreeze and VitriThaw (VitriFreeze/Thaw) are a set of ready-to-use media for vitrification and thawing of human embryos.

*For professional use only.*

**BACKGROUND**

Vitrification which is preservation at extremely low temperatures without freezing can be more favourable than slow cooling (Kuleshova and Lopata, 2002). Due to variable results after application of slow-freezing methods for blastocysts, vitrification was introduced as an alternative approach.

The success rates of vitrification have been increased with the use of ultra-rapid vitrification procedures. Since several years, vitrification of blastocysts and embryos using different embryo carriers has resulted in many pregnancies.

In order to accomplish aseptic vitrification of blastocysts, the HSV High security vitrification kit (Cryo Bio System; available from FertiPro N.V.) was developed (Vanderzwalmen et al, 2000). The tip of the HSV straw is designed to hold the blastocysts in a very small volume of vitrification solution allowing fast cooling and thawing rates as compared to those achieved by immersion of sealed 0.25ml plastic straws into liquid nitrogen. An alternative, aseptic method is now also available in the shape of the VitriSafe straw.

**COMPOSITION**

VitriFreeze/Thaw are DMSO/ethylene glycol based vitrification media that also contain PBS, sucrose, Ficoll and human serum albumin (10-20g/Liter). VitriFreeze/Thaw do not contain antibiotics.

**MATERIAL INCLUDED WITH THE KIT**

**VitriFreeze™ kit (VF\_KIT1)**

- » 1 vial of VitriFreeze Pre-incubation medium (5ml)
- » 1 vial of VitriFreeze Freezing medium 1 (1ml)
- » 1 vial of VitriFreeze Freezing medium 2 (1ml)

**VitriThaw™ kit (VT\_KIT1)**

- » 1 vial of VitriThaw Thawing medium 1 (5ml)
- » 1 vial of VitriThaw Thawing medium 2 (1ml)
- » 1 vial of VitriThaw Thawing medium 3 (1ml)
- » 1 vial of VitriThaw Thawing medium 4 (1ml)

The media should be used in the order displayed above (the bottles may be in a different order in the kit) and can be used for approximately 4 procedures.

**MATERIAL NOT INCLUDED WITH THE KIT**

- » well dishes (e.g. Nunc 144 444)
- » freezing tank with liquid nitrogen
- » water bath (able to hold 37°C)
- » attenuated pipettes
- » forceps
- » vitrification device (HSV device, Cryo Bio System)
- » LAF-bench (ISO Class 5)
- » microscope
- » lab timer

**VITRIFREEZE/THAW AND EMBRYOCULTURE**

VitriFreeze/Thaw can be used in combination with FertiCult media (Flushing, IVF and GAIN medium) before freezing and after thawing.

**PRODUCT SPECIFICATIONS**

- » Chemical composition
- » pH between 7,20 – 7,40
- » Osmolality:
  - › VitriFreeze Pre-incubation medium: 270-290mOsm/kg
  - › VitriThaw Thawing medium 4: 270-290mOsm/kg
  - › VitriThaw Thawing medium 1: 805-850 mOsm/kg
  - › VitriThaw Thawing medium 2: 535-565 mOsm/kg
  - › VitriThaw Thawing medium 3: 405-435 mOsm/kg
- » Sterility: SAL 10<sup>-3</sup>
- » Endotoxins: < 0,25 EU/ml
- » Mouse Embryo Assay (blastocysts after 96h) ≥ 80%
- » Use of Ph Eur or USP grade products if applicable
- » The certificate of analysis and MSDS are available upon request

**PRE-USE CHECKS**

- » Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- » Do not use the product if seal of the container is opened or defect when the product is delivered.

**STORAGE INSTRUCTIONS**

- » Store between 2-8°C.
- » Do not freeze before use.
- » Keep away from (sun)light.
- » The products can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- » Do not use after expiry date.
- » Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).

**WARNINGS AND PRECAUTIONS**

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination.

Only for the intended use. The long term safety of embryo vitrification on children born following this procedure is unknown.

**METHOD**

Ensure all media are well mixed before use. We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

**Preliminary steps**

In a 4-well dish fill the first well with 300 µl of Pre-incubation medium, the second with VitriFreeze 1 (250 µl) and the third with VitriFreeze 2 solution (250 µl). Next open as many packs of HSV devices as will be required for the vitrification step, taking into account that 1 HSV device can hold up to 2 embryos. Conveniently place the separate parts of the HSV device on the workbench for easy access later in the procedure. Up to 5 vitrification cycles (of the same patient!) can be performed with one media set-up. Do not use the same media for different patients!

**Freezing preparation**

Transfer the embryos from the blastocyst cell culture medium in to each of the VitriFreeze solutions using the following scheme:

Stage	Pre-incub.	Vitri Freeze 1	Vitri Freeze 2	Temperature
Early Blastocyst Morulae	2'	2'	30"	Room temperature
Blastocyst – expanded Blastocyst	2'	3'	30"	37°C
Blastocyst – expanded Blastocyst + artificial shrinkage*	2'	2'	30"	Room temperature

\* Before starting the vitrification procedure, in order to reduce the negative effect of the blastocoeal, expanded blastocysts should be collapsed by reducing artificially with a glass pipette the volume of the blastocoeal (Vanderzwalmen et al, 2002, Son et al., 2003 Hiraoka 2004).

**Vitrification**

- Using an attenuated pipette or an equally suitable device, deposit maximum 2 blastocysts in a volume of approximately 0,3µl of VitriFreeze 2, in the gutter of the tip of the vitrification straw.
- Place the vitrification straw in the outer sheath and seal it as indicated in the instructions for use of the HSV device.
- Plunge the sealed device into the liquid nitrogen.

**Thawing**

- In a 4-well dish, fill the first well with 1 ml VitriThaw Thawing medium 1, the second with Thawing 2 (250 - 300 µl), the third with Thawing 3 (250 - 300 µl) and the fourth with Thawing 4 (250 - 300 µl). Pre-heat the media at 37°C.
- Remove the vitrification straw from the outer sheath as indicated in the instructions for use of the HSV device.
- Immediately plunge the vitrification straw into pre-heated VitriThaw Thawing medium 1 (37°C) and leave in thawing 1 for 3 minutes.
- Transfer into VitriThaw Thawing medium 2 (37°C) and leave in this medium for 2 minutes.
- Transfer into VitriThaw Thawing medium 3 (37°C) and leave in this medium for 2 minutes.
- Finally transfer into VitriThaw Thawing medium 4 (37°C) and wash for at least 1 minute.
- Transfer into blastocyst culture medium for continued cell culture.

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