VitriFreeze™

VitriThaw™

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BIBLIOGRAFIE

1. Hiraoka K., Kinutami M et al. (2004). Blastocoele collapse by micropipetting prior to vitrification gives excellent survival and pregnancy outcomes for human day 5 and 6 expanded blastocysts. Hum. Reprod., 19, 2884–2888.
2. Kuleshova, L. and Lopata, A. (2002) Vitrification can be more favourable than slow cooling. Fertil. Steril., 78, 449-454.
3. Son, W.Y., Yoon, S.H., Yoon, H.J., Lee, S.M. and Lim, J.H. (2003) Pregnancy outcome following transfer of human blastocyst vitrified on electron microscopy grids after induced collapse of the blastocoel, Hum. Reprod., 18, 137-139.
4. Vanderzwalmen, P., Bertin, G., Debauche, C.H., Standaert, V. And Schoysman, R. (2000) In vitro survival of metaphase II oocytes and blastocysts after vitrification in an Hemi Straw (HS) system, Fertil. Steril., 74 (Suppl.), 215.
5. Vanderzwalmen, P., Bertin, G., Debauche, C.H., Standaert, V., van Roosendaal, E., Vandervorst, M., Bollen, N., Zech, H., Mukaida, T., Takahashi, K. and Schoysman, R. (2002) Births after vitrification at morula and blastocyst stages: effects of artificial reduction of the blastocoelic cavity before vitrification. Hum. Reprod., 17, 744-751.
6. Vanderzwalmen, P., Bertin, G., Debauche, C.H., Standaert, V., Bollen, N., van Roosendaal, E., Vandervorst, M., Schoysman, R. and Zech, H. (2003) Vitrification of human blastocyst with the Hemi-Straw carrier: application of assisted hatching after thawing. Hum. Reprod., 7, 1-8.
7. S. Shen, X. Yang, M.J. Abeyta, S.L. Benedict, M.P. Rosen, M.I. Cedars (2008). Vitrification using a closed cryo-top system significantly improves day 3 embryo survival compared to slow freezing. Fertil. Steril., 90 (Suppl 1.), S293-S293.
8. Zech, N.H., Lejeune, B., Zech, H., Vanderzwalmen, P. (2005). Vitrification of hatching and hatched human blastocysts: effect of an opening in the zona pellucida before vitrification. Reproductive BioMedicine Online, 11 (No 3), 355-36.
9. Jelinkova, L., Ditzel, N., Reeka, N., Gagsteiger, F., Moosova, M., Pavelkova, J., Rezabek, K. (2006). The removal of ZP from blastocysts before vitrification increased their survival and also their implantation rates. Abstracts of the 22nd Annual Meeting of the ESHRE, Prague, Czech Republic, 18–21 June 2006.
10. Shebl, O., Ebner, T., Sommergruber, M., Sir, A., Tews, G. (2009). Cryopreserved blastocysts have a lower implantation but an equal live birth rate as compared to fresh blastocysts of the same quality – a case – control study. Acta Obstetrica et Gynecologica Scandinavica 88 (No. 8), 944-947.
11. Ebner, T., Vanderzwalmen, P., Shebl, O., Urdl, W., Moser, M., Zech, N.H., Tews, G. (2009). Morphology of vitrified/warmed day-5 embryos predicts rates of implantation, pregnancy and live birth. Reproductive BioMedicine Online, 19 (No. 1), 72-78.

TECHNISCHE ONDERSTEUNING

FertiPro N.V.
Industriepark Noord 32
8730 Beernem, Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

CE
0344