VitriFreeze ES™ VitriThaw ES™

VitriFreeze ES™ VitriThaw ES™

VitriFreeze ES™ VitriThaw ES™

Universal media for vitrification and warming of human embryos: zygote, cleavage stage, blastocyst

STERILE A

Sterilised by sterile filtration. Doc. reference: FP09 I46 02 R01 D.7 Update: 04.03.2019

INTENDED USE

VitriFreeze ES™ and VitriThaw ES™ (VitriFreeze/Thaw ES™) are a set of ready-to-use media for vitrification and warming of human embryos. The Early stage vitrification kit and protocol are meant to be used with human embryos between 2PN and blastocyst stage.

For professional use only.

BACKGROUND

Vitrification is a cryopreservation procedure where a liquid solution is converted into an amorphous solid, free of any crystalline structure (Rall and Fahy 1985). This technique can be more favourable than slow cooling (Stehlik 2005). Ultra-rapid vitrification of zygotes and embryos using 'open' carrier devices such as the 'Hemi-Straw' or 'VitriPlug', which allow direct contact with liquid nitrogen, has resulted in many births of healthy babies (Vanderzwalmen, 2003). Due to European regulations defining medical safety requirements for cryopreservation of human cells, hermetically closed (aseptic) containers were developed which avoid direct contact between the embryo and the liquid nitrogen during cooling and long term storage. For this purpose the HSV (High Security Vitrification kit - Cryo Bio System) and the VitriSafe plug (MTG) were developed (Vanderzwalmen, 2009). Both devices consist of an inner straw that contains a gutter in which a small volume of cryoprotectant solution, with one or two embryos, is deposited. This inner straw is inserted into an outer, protective straw that is sealed before immersion into liquid nitrogen. Due to thermo isolation however, the cooling rate in these devices is reduced compared to that in 'open' carrier devices. Therefore more cryoprotectant has to penetrate the cells in order to guarantee an intracellular state. The VitriFreeze ES™ kit is designed to allow sufficient amounts of cryoprotectant to enter the embryo at different stages of development. During warming, VitriThaw ES™ kit is designed to gradually remove cryoprotectants.

COMPOSITION

VitriFreeze/Thaw ES™ are DMSO/Ethylene Glycol based vitrification media that also contain PBS, sucrose, Ficoll and human serum albumin (10-20g/liter). VitriFreeze/Thaw ES™ media do not contain antibiotics.

MATERIAL INCLUDED WITH THE KIT

- VitriFreeze ES™ kit (VF_KIT1_ES) contains one bottle each of the following media:
 - 5ml VitriFreeze ES™ - Pre-incubation medium ("VPI")
 - 1ml VitriFreeze ES™ - Freezing medium 1 (5% DMSO - 5% EG) ("VF1")
 - 1ml VitriFreeze ES™ - Freezing medium 2 (10% DMSO - 10% EG) ("VF2")
 - 1ml VitriFreeze ES™ - Freezing medium 3 (20% DMSO - 20% EG) ("VF3")

- VitriThaw ES™ kit (VT_KIT1_ES) contains one bottle each of the following media:
 - 5ml VitriThaw ES™ - Thawing medium 1 ("VT1")
 - 3.2ml VitriThaw ES™ - Thawing medium 2 ("VT2")
 - 1ml VitriThaw ES™ - Thawing medium 3 ("VT3")
 - 1ml VitriThaw ES™ - Thawing medium 4 ("VT4")
 - 1ml VitriThaw ES™ - Thawing medium 5 ("VT5")

The media should be used in the order displayed above (the bottles may be in a different order in the box) and can be used for approximately 3-4 procedures.

MATERIAL NOT INCLUDED WITH THE KIT

- Well dishes (e.g. Nunc 144 444)
- Freezing tank with liquid nitrogen
- Water bath (able to hold 37°C)
- Attenuated pipettes
- Forceps
- Vitrification device (HSV device, VitriSafe)
- LAF-bench (ISO Class 5), Microscope, Lab timer

VITRIFREEZE/THAW ES AND EMBRYOCULTURE

VitriFreeze/Thaw ES™ can be used in combination with GAIN™ media/FertiCult™ media (Flushing, IVF) before freezing and after thawing.

PRODUCT SPECIFICATIONS

- Chemical composition
- pH between 7,20 – 7,40
- Osmolality:
 - VPI: 270 – 295 mOsm/kg (release criteria: 270 – 290)
 - VT2: 805 – 865 mOsm/kg (release criteria: 805 – 850)
 - VT3: 535 – 565 mOsm/kg
 - VT4: 405 – 435 mOsm/kg
 - VT5: 270 – 295 mOsm/kg (release criteria: 270 – 290)
- Sterility: SAL 10⁻³
- Endotoxins: < 0,25 EU/ml
- Mouse Embryo Assay (blastocysts after 96h) ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- The certificate of analysis and MSDS are available upon request

PRE-USE CHECKS

- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defect when the product is delivered.
- VitriFreeze ES™ - Freezing medium 2 may contain small salt precipitates which do not have an impact on product performance/safety.

STORAGE INSTRUCTIONS

- Store between 2-8°C.
- Do not freeze before use.
- Keep away from (sun)light.
- The products can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Do not use after expiry date.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).

WARNINGS AND PRECAUTIONS

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination. Only for the intended use. The long term safety of embryo vitrification on children born following this procedure is unknown.

METHOD

Ensure all media are well mixed before use. We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps

- Up to 5 vitrification cycles (of the same patient) can be performed with one media set-up prepared as indicated below. Do not use the same media for different patients!
- Open as many packs of vitrification devices as will be required for the vitrification step, taking into account the 1 device can hold 1 or 2 embryos (check the instructions for the device you are using). Conveniently place the separate parts of the device on the workbench for easy access later in the procedure.
- Cooling procedure: In a 4-well dish fill:
 - 250-300µl VPI
 - 250-300µl VF1
 - 250-300µl VF2
 - 250-300µl VF3
- Warming procedure: in a 6-well dish, fill:
 - 500-800µl VT1 (with closed device)
 - 250-300µl 1:1 dilution of VT1 and VT2 (with closed device)
 - 250-300µl VT2 (with closed device)
 - 500-800µl VT2 (with open device)
 - 250-300µl VT3
 - 250-300µl VT4
 - 250-300µl VT5

Cooling protocol using a closed (aseptic) device

Warm all media of the kit to room temperature (20-25°C) before use.

Development stage	VPI	VF1	VF2	VF3
Zygote-2-cell	2'	5'-10'	4'-5'	40"-60"
4-8 cell	2'	5'-7'	4'	40"-60"
Morula	2'	5'-7'	4'	40"-60"
Early Blastocyst	2'	5'-7'	4'	40"-60"
Expanded blastocyst	2'	5'-10'	4'	40"-60"

Note: The complete process of placing the embryo in VitriFreeze ES™ Freezing medium 3, loading the embryo on the vitrification device, inserting the device in the outer straw (if applicable) and sealing should not take longer than 60 seconds before plunging the device into the liquid nitrogen. In case the process has taken longer than 60 seconds, make a note of this to analyze the effect on the results afterwards.

Cooling protocol using an open device

Warm all media of the kit to room temperature (20-25°C) before use.

Development stage	VPI	VF1	VF2	VF3
Zygote-2-cell	2'	2'	3'	30"-40"
4-8 cell	2'	2'	3'	30"-40"
Morula	2'	2'	3'	30"-40"
Early Blastocyst	2'	2'	3'	30"-40"
Expanded blastocyst*	2'	2'	4'	30"-40"

* If artificial shrinkage is used (Vanderzwalmen, 2002; Son,2003)

Note: The complete process of placing the embryo in VitriFreeze™ ES Freezing medium 3, loading the embryo on the vitrification device, inserting the device in the outer straw (if applicable) and sealing should not take longer than 60 seconds before plunging the device into the liquid nitrogen. In case the process has taken longer than 60 seconds, make a note of this to analyze the effect on the results afterwards.

Warming protocol using a closed device

Warm all media of the kit to room temperature (20-25°C) before use. (Alternatively and as indicated in previous version of the instructions for use, VitriThaw ES™ - Thawing medium 1 can be warmed to 37°C.)

Development stage	VT1	VT1/2*	VT2	VT3	VT4	VT5
Zygote-2-cell	1'	1'	1'-2'	2'-4'	2'-4'	Wash for 1'-2' before transfer to culture medium
4-8 cell	1'-1'30"	1'-1'30"	1'-2'	2'-4'	2'-4'	
Morula	1'-1'30"	1'-1'30"	1'-2'	2'-4'	2'-4'	
Early Blastocyst	1'-1'30"	1'-1'30"	1'-2'	2'-4'	2'-4'	
Expanded blastocyst	1'-1'30"	1'-1'30"	1'-2'	2'-4'	2'-4'	

* Mix 1 part VT1 with 1 part VT2

Warming protocol using an open device

Warm all media of the kit to room temperature (20-25°C) before use. (Alternatively and as indicated in previous version of the instructions for use, VitriThaw ES™ - Thawing medium 2 can be warmed to 37°C.)

Development stage	VT2	VT3	VT4	VT5
Zygote-2-cell	2'	2'-4'	2'-4'	Wash for 1'-2' before transfer to culture medium
4-8 cell	2'	2'-4'	2'-4'	
Morula	2'	2'-4'	2'-4'	
Early Blastocyst	2'	2'-4'	2'-4'	
Expanded blastocyst	2'	2'-4'	2'-4'	

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